

DECEMBER 2016

## Medical Research Symposium & Pacific Partnership in conjunction with

## Kuantan Research Day 2016

International Research Network: The Way Forward

## **ABSTRACTS**

Organized by:



Kulliyyah of Medicine, IIUM Research Management Centre, Kuantan









Basic Health Sciences Oral

Abstract ID: 34

## Epigenetic methylation status of P16, MGMT and SPOCK2 in diffuse Large B cell lymphoma

**Lailatul Jalilah**<sup>a</sup> | Aziah Ismail<sup>b</sup> | Norlelawati A. Talib<sup>c</sup> | Naznin Muhammad<sup>c</sup> | Faezahtul Arbaeyah Hussain<sup>d</sup> | Norafiza Zainuddin<sup>a</sup>

<sup>a</sup>Department of Biomedical Science, Kulliyyah of Allied Health Sciences, International Islamic University Malaysia

bInstitute for Research in Molecular Medicine (INFORMM), Universiti Sains Malaysia

<sup>c</sup>Department of Pathology and Laboratory Medicine, Kulliyyah of Medicine, International Islamic University Malaysia

<sup>d</sup>Department of Pathology, School of Medical Sciences, Universiti Sains Malaysia

Introduction: Epigenetic methylation has been implicated in the pathogenesis of diffuse large B cell lymphoma (DLBCL). This study investigated the methylation status of p16, MGMT and SPOCK2. Aberrantly methylated p16 and MGMT have been linked to DLBCL, but not SPOCK2. p16 inhibits cyclin-dependent kinase, which results in retinoblastoma phosphorylation and blockage of cell cycle at G1 phase. MGMT removes alkyl adduct at O6guanine, thus preventing lethal cross-links. SPOCK2, an extracellular chondroitin and heparin sulfate proteoglycans, abolishes the inhibition of membrane-type 1-matrix metalloproteinase which might enhance the angiogenesis. The absence of SPOCK2 methylation was therefore hypothesized in the majority of cases in this study. Methods: Extracted DNA from 88 formalin-fixed paraffin-embedded (FFPE) tissues of DLBCL were subjected to bisulfite conversion followed by methylation-specific PCR (MSP) analysis for p16, MGMT and SPOCK2 methylation. p16 methylation was also quantified in 16 samples through pyrosequencing assay. Results: p16 methylation was observed in 65/88 (74%) samples by MSP. Pyrosequencing detected p16 methylation in all 16 samples ranging from 18% to 81%. MGMT methylation was detected in all 88 (100%) cases. Methylated SPOCK2 was found in 83 (94.3%) samples. There was a significant association between p16 methylation status with patients above 50 years of age (p=0.04). Conclusions: These preliminary discoveries may serve as a good platform in order to gain a comprehensive overview on the epigenetics contribution in the pathogenesis of DLBCL. Pyrosequencing is a robust tool in detecting and quantifying methylation.

KEYWORDS: DLBCL, epigenetics, MSP, pyrosequencing